

# Chromatography

## QuEChERS Bulk Sorbents

GVS provides superior quality QuEChERS bulk sorbents which have been verified by our lab. You can choose the ratio according to the needs of the experiment.

### Features

- Supply 2 mL or 15 mL clean-up tubes
- Suitable for AOAC 2007, EN 15662 standards, etc.



### Ordering information

Product Code	Sorbent	Specification	Qty.
SPESORB00PSA2100A	PSA	Carbon Content: 8%, Surface area: 480 m <sup>2</sup> /g, Particle size: 50-75 µm, Pore size: 70 Å	100 g
SPESORB00C181100A	C18	Carbon Content: 17.6%, Surface area: 300 m <sup>2</sup> /g, Particle size: 40-75 µm, Pore size: 70 Å	100 g
SPESORB00GCB1050A	Carb-GCB	Surface area: 100 m <sup>2</sup> /g, Particle size: 100-300 mesh	50 g

## QuEChERS Ceramic Homogenizers

GVS Ceramic Homogenizers can be used in QuEChERS extraction kit and clean-up kit to improve extraction recovery and reproducibility.

### Features

- Inert ceramic material, no impurities dissolution
- Shorten sample extraction time and reduce labor cost
- Improve recovery rate and reproducibility of sample extraction



### Ordering information

Product Code	Description	Qty.
QUEBCER02YYA	Ceramic Homogenizers, 50 mL	100 Pcs/Box
QUEBCER15YYA	Ceramic Homogenizers, 15 mL	100 Pcs/Box
QUEBCER50YYA	Ceramic Homogenizers, 2 mL	200 Pcs/Box