


## Cell Culture Flask



- According to different surface treatment processes, cell culture flasks can be divided into hydrophilic cell culture flask and hydrophobic cell culture flask.
- Under the same environment, the clone formation rate and growth generation of cells cultured in the cell culture flasks have reached the index of similar international products.
- High quality polystyrene (PS) materials, USP VI medical grade.
- Plasma surface treatment can maximize adhesion of cells. Culturing for 48h, Vero cell growing condition is good as international famous culture flask.
- Innovative angled neck design provides good pipette and cell scraper access.
- The whole production is in the C-class purification

workshop. According to ISO13485 quality system management and the use of automatic equipment production, no direct contact with personnel, good product consistency, and small differences between batches.

- Sterile, DNase-free, RNase-free, endotoxin content < 0.1 EU/mL.
- According to different culture requirements, there are vent cap and seal cap optional.
- Upgrade the flask angle and bottom arc design to easily obtain a complete growth surface and make subsequent operations easier.



## Cell Culture Flask

### Ordering information

TC Treated Flask with Seal Cap					
Product Code	Size	Culture area [cm <sup>2</sup> ]	Cap Type	pcs/pack	pcs/case
CELCUFG006025S	25 cm <sup>2</sup>	25	Seal cap	10	200
CELCUFG006075S	75 cm <sup>2</sup>	75	Seal cap	5	100
CELCUFG006175S	175 cm <sup>2</sup>	175	Seal cap	5	40
CELCUFG006225S	225 cm <sup>2</sup>	225	Seal cap	5	25

TC Treated Flask with Vent Cap					
Product Code	Size	Culture area [cm <sup>2</sup> ]	Cap Type	pcs/pack	pcs/case
CELCUFG004025S	25 cm <sup>2</sup>	25	Vent cap	10	200
CELCUFG004075S	75 cm <sup>2</sup>	75	Vent cap	5	100
CELCUFG004175S	175 cm <sup>2</sup>	175	Vent cap	5	40
CELCUFG004225S	225 cm <sup>2</sup>	225	Vent cap	5	25

Non-TC Treated Flask with Seal Cap					
Product Code	Size	Culture area [cm <sup>2</sup> ]	Cap Type	pcs/pack	pcs/case
CELCUFG007025S	25 cm <sup>2</sup>	-	Seal cap	10	200
CELCUFG007075S	75 cm <sup>2</sup>	-	Seal cap	5	100
CELCUFG007175S	175 cm <sup>2</sup>	-	Seal cap	5	40
CELCUFG007225S	225 cm <sup>2</sup>	-	Seal cap	5	25

Non-TC Treated Flask with Vent Cap					
Product Code	Size	Culture area [cm <sup>2</sup> ]	Cap Type	pcs/pack	pcs/case
CELCUFG005025S	25 cm <sup>2</sup>	-	Vent cap	10	200
CELCUFG005075S	75 cm <sup>2</sup>	-	Vent cap	5	100
CELCUFG005175S	175 cm <sup>2</sup>	-	Vent cap	5	40
CELCUFG005225S	225 cm <sup>2</sup>	-	Vent cap	5	25

